

Journal Club – Optics

11/05/2020

SCIENCE ADVANCES | RESEARCH ARTICLE

BIOPHYSICS

Ultraprecise single-molecule localization microscopy enables in situ distance measurements in intact cells

Simao Coelho^{1*}, Jongho Baek^{1*}, Matthew S. Graus¹, James M. Halstead¹, Philip R. Nicovich², Kristen Feher¹, Hetvi Gandhi¹, J. Justin Gooding³, Katharina Gaus^{1†}

<https://advances.sciencemag.org/content/6/16/eaay8271>

Scope/Claim of the paper



- direct distance measurements between molecules in intact cells on the scale between 1 and 20 nm
 - biological relevant length scales
- stabilization of <1 nm and localization precision of ~ 1 nm
 - “drift is reduced below the photon-limited localization precision and causes the fluorescent positions to be prealigned”

Intro


- engineering solution for SMLM, termed Feedback SMLM
 - SMLM: long duration of experiment (6h), but also short binding events (200ms)
- No post-acquisition processing including drift correction, grouping, filtering, averaging, or summation
 - Post-acquisition processing has problems with high densities of detections such as proteins in clusters with and with rare events

Method

- real-time drift correction in 3D, using three types of corrections
 - nonfluorescent fiducials are used outside the field of view of the camera for SMLM
 - autonomous optical feedback loop in the emission path
 - Reduction of the variation registered across the EMCCD







Setup



Setup

EMCCD with $40\mu\text{m} \times 40\mu\text{m}$ FoV (pixel size, $80 \times 80 \text{ nm}$)
100 \times 1.49 NA TIRF objective (Nikon, CFI Apochromat)



Fiducials


- 3 μm polystyrene beads that are illuminated with an infrared LED to create diffraction rings
 - Provide information of both the x-y and z position
 - Look-up-table for the z direction
- CMOS camera (112 μm \times 70 μm FoV, about fivefold larger than the fluorescence camera)
- stage corrections at 15 Hz, independent of imaging
- stabilization of 0.4 and 1 nm (SD) in the lateral and axial directions

Feedback loop in the emission

- White LED, integrated in the microscope body, creates an optical fiducial on the EMCCD
 - position of LED is localized with a precision of 0.05 nm
 - piezoelectric mirror is used to account for long-term drift of the microscope (i.e., mechanical instabilities)
 - reduces the image drift to 0.22 nm
 - LED spot was positioned on the periphery of the EMCCD

EMCCD characterization and chromatic correction

- “nanohole array”
 - holes with a diameter of 100 nm nanofabricated on an aluminum-coated coverslip
 - filled with red (Alexa Fluor 647) and green (Alexa Fluor 488) dyes at $\sim 100 \mu\text{M}$ concentration



Only chromatic aberration or also other aberrations corrected?


Additional info

- environmental control box was installed around the optical setup to control temperature fluctuations (peltier elements)
 - custom-built
 - temperature and humidity have a standard deviation of 0.02°C and 0.88% , respectively
- data acquisition and instrument control were performed using an Xbox One controller

DNA-PAINT

- No post-acquisition processing
- ~1 nm localization precision
- No 3D localization, statistics?


DNA origami w/ 3 docking sites
 “Standard” is actually a commercial Zeiss microscope



actin filaments


- identify single filaments and determine their widths as 5 to 9nm (again no post-acquisition processing)
- No residual drift when using gold nanoparticles

phalloidin in a COS-7 cell



active and resting T cells in situ

- DNA-PAINT acquisitions of CD45 and phosphorylated CD3 ζ
- Difference between dist in C and D



Conclusion

- Claims could be better backed by data
 - Only 1 image of the DNA origami. A bit the same for all other analysis/samples.
 - Distance measurements are confuse. The used system is really not ideal.
- How does the LED spot on the EMCCD look like?
- Not easy to mount an LED in a normal microscope body